

# Keysight Technologies

## In Vitro Complex Shear Modulus of Bovine Muscle Tissue (steak)


### Application Note

#### Abstract

Dynamic instrumented indentation provides a way to measure the mechanical properties of soft biological tissue that is relevant, accurate, and easy. In this work, we measured the complex shear modulus of bovine muscle tissue, submerged in saline at bovine body temperature (38°C). Along the grain, the shear modulus was  $24.4 \pm 13.9$  kPa. Orthogonal to the grain, the shear modulus was  $11.4 \pm 2.9$  kPa. The loss factor was highly consistent and independent of testing direction:  $\tan \delta = 0.34 \pm 0.035$ .

#### Introduction

Being able to measure the mechanical properties is crucial to understanding soft living tissue. First, there is a reciprocal relationship between mechanical properties and function. That is, properties affect function, but function can also affect properties by means of adaptation. Thus, knowledge of mechanical properties leads directly to knowledge of function. Further, gel is often used as a substitute for living tissue [1]. In order to tailor the gel to mimic the tissue, one must first know the mechanical properties of the tissue. Thus, the purpose of this work was to demonstrate the use of instrumented indentation to measure the complex shear modulus of soft tissue under physiologically relevant



conditions—submerged in saline at body temperature.

Figure 1. Schematic of skeletal muscle tissue (National Cancer Institute, public domain).

We chose bovine muscle tissue (steak) as an exemplary material. Figure 1 shows a schematic of skeletal muscle tissue. The “grain” of the steak runs in the long direction of the muscle tissue; this is the direction in which the muscle contracts and relaxes. A cut the grain of the steak reveals the same cross-section shown schematically in Figure 1. The proper term for a “grain” is the , which is a bundle of individual muscle cells or fibers. Given this structure of muscle tissue, we should expect the properties of the tissue to be anisotropic.

Dynamic instrumented indentation returns the stiffness ( $S$ ) and damping ( $C\omega$ ) of the contact at the prescribed frequency. Using Sneddon’s elastic contact theory<sup>2</sup> as later developed by Oliver, Pharr, and Brotzen<sup>3</sup> and Herbert et al.<sup>4</sup>,

the real and imaginary components of the complex shear modulus ( $G^* = G' + iG''$ ) are calculated as

$$G' = S(1-\nu)/(2D), \text{ and}$$

$$G'' = C\omega(1-\nu)/(2D),$$

where  $\nu$  is the Poisson's ratio of the material (assumed 0.5 for soft tissue) and  $D$  is the contact diameter (i.e. the diameter of the punch face). Traditionally, the shear loss modulus ( $G''$ ) is not reported as an absolute value, but in relation to the shear storage modulus. The loss factor, which characterizes the damping relative to stiffness, is calculated directly as

$$\tan\delta = G''/G' = C\omega/S.$$

## Experimental Method

The Keysight Technologies, Inc. G200 NanoIndenter was used for all testing. The system was configured with a DCM II actuator, the CSM option, and the hot-stage option. The CSM option allowed the superposition of an oscillating force, and the hot-stage option allowed the samples to be tested at body temperature. The DCM II actuator was fitted with a flat-ended cylindrical punch ( $D=101.1\mu\text{m}$ ) in order to generate a known and constant contact.




Figure 2. Steak samples as prepared for testing (a) cut across the grain, and (b) cut with the grain.

Two samples of steak were tested in this work: one cut across the grain (Figure 2a) and one cut with the grain (Figure 2b). Each kind of sample was prepared in the following way. The night before testing, the tissue was rough-cut and frozen. A disk ( $\approx 8\text{mm}$ ,  $\approx 8\text{mm}$ ) was cut from the frozen material and quickly adhered to the base of the sample well<sup>1</sup> using 5-minute epoxy. Once the epoxy had set, the well was filled with saline, and the sample holder was attached to the hot stage. Finally, a

microscope slide was gently laid over the sample to keep it from drying while the sample equilibrated at the testing temperature ( $38^\circ\text{C}$ ). The equilibration time was about 30 minutes. Figure 3 shows a sample ready for testing.

It should be noted that the testing direction is orthogonal to the plane of the cut. Thus, the sample which was cut across the grain (Figure 2a) was actually tested in the direction of the grain—that is, in the direction of the muscle fibers illustrated in Figure 1. The sample which was cut with the grain was actually tested perpendicular to the grain. Results are reported according to the direction of the *test*, not the cut.

The NanoSuite test method “G-Series DCM CSM Flat Punch Complex Modulus, Gel” was used for all testing. Each test with this method comprised the following steps:

1. Approach the surface, oscillating the indenter at its natural frequency (110Hz), until contact is detected,
2. Switch to the user-prescribed testing frequency (10Hz),
3. Apply pre-test compression ( $7\mu\text{m}$ ),
4. Sense tissue stiffness and damping, and
5. Withdraw the indenter and move to the next test site.



Figure 3. Steak sample, submerged and ready for testing at body temperature. Microscope slide was used to seal the sample and keep it moist until just prior to testing.

1. In order to allow the samples to be submerged during testing, we used washers to build a “well” on the sample base plate for the hot-stage. Two washers, built up from the base plate with JB Weld SteelStik, provided sufficient depth (8mm).

With this method, twenty-five different sites were tested on each sample; individual test sites were separated by at least  $400\mu\text{m}$ . Tests were performed in red tissue only. Each test took about two minutes. In order to keep the samples from drying out, testing was paused after every 5 tests in order to add 1–2 drops of saline. Others have solved the problem of drying by integrating an intravenous drip [5].

## Results and Discussion

Testing a sample submerged in saline poses no problems. Although the system indeed senses contact with the saline, the signature is different from that of contact with tissue (Figure 4). The attraction between the tip and the fluid causes a momentary decrease in the stiffness of the system. Thus, the interaction between the tip and the fluid manifests as a slight, but detectable, increase in the phase angle. Subsequent contact with the sample causes an increase in the stiffness of the system, thus manifesting as a sharp decrease in the phase angle. This extreme sensitivity in the phase angle is attained by oscillating the indenter at the natural frequency of the instrument (110Hz) during the approach.

Figure 5 shows the results for shear modulus. The muscle tissue is significantly stiffer in the direction of the grain ( $G' = 24.4 \pm 13.9 \text{ kPa}$ ) than perpendicular to




Figure 4. Signature of contact with saline and tissue.

it ( $G' = 11.4 \pm 2.9 \text{ kPa}$ ). This is not surprising since the muscle tissue naturally acts in the direction of the grain to alternately exert and relax force.

The point-to-point variation in modulus is high, especially in the direction of the grain, but this is to be expected for biological materials. The observed variation is due to true point-to-point variation in properties, not measurement uncertainty<sup>2</sup>. Indeed, when we test uniform gels of comparable moduli, the point-to-point variation is quite low [6]. In the direction of the grain, this high degree of variation is likely due to the complex structure exposed by the cut across the grain (as illustrated by

the cross-section of Figure 1). Perpendicular to the grain, the variation is smaller. Again, this is not surprising, because the structure is less complex in this direction.

The uniformity and consistency of the loss factor is quite surprising, however (Figure 6). The ability of the tissue to damp energy, as quantified by the loss factor is rather independent of the direction or location ( $\tan \delta = 0.34 \pm 0.035$  for both directions). This surprising result begs for a phenomenological explanation. The value of the loss factor (0.34) means that at 10Hz, the capacity of the tissue to damp out energy is one third of its capacity to store energy elastically. This is physiologically realis-




Figure 5. Shear modulus of bovine muscle tissue at  $38^\circ\text{C}$ , 10Hz (N = 25).




Figure 6. Loss factor of bovine muscle tissue at  $38^\circ\text{C}$ , 10Hz (N = 25).

2. At 10Hz, the uncertainty in the stiffnesses measured by the DCM II is less than  $0.1\text{N/m}$ . In this work, the contact stiffnesses were  $5\text{--}10\text{N/m}$ . Thus, we expect 1–2% variation due to measurement uncertainty.

tic and only measurable by keeping the sample submerged during testing.

## Conclusions

The values for complex shear modulus measured in this work confirm the anisotropic nature of muscle tissue and provide a reliable basis for relating other comparable measurements. These extraordinary measurements are made possible by unique hardware and software features of the Keysight G200 NanoIndenter. The ability to oscillate the indenter makes the system sensitive to the slightest changes, and the software automatically responds

to such changes to direct the progress of the experiment according to the experimenter's prescription. We anticipate that these techniques will be used by others to illuminate fundamental relationships between structure, history, and mechanical properties in biological tissues.

## Acknowledgement

This work was funded in part by the Bulgarian Ministry of Education, Youth and Science, Human Resources Development Operational Program, Project BG051PO001/3.3-05-0001 "Science and Business".

## References

1. Jain, S.K., Bhattacharayya, C.N., Badonia, B., and Singh, R.P., "Study of Unusual Phenomenon of Contact Firing on Gelatine Block Using .38 Special Revolver—Forensic Importance," *Forensic Science International* 133(3), 183–189, 2003.
2. Sneddon, I.N., "The Relation Between Load and Penetration in the Axisymmetric Boussinesq Problem for a Punch of Arbitrary Profile," *Int. J. Eng. Sci.* 3(1), 47–57, 1965.
3. Pharr, G.M., Oliver, W.C., and Brotzen, F.R., "On the Generality of the Relationship among Contact Stiffness, Contact Area, and Elastic-Modulus during Indentation," *Journal of Materials Research* 7(3), 613–617, 1992.
4. Herbert, E.G., Oliver, W.C., and Pharr, G.M., "Nanoindentation and the Dynamic Characterization of Viscoelastic Solids," *Journal of Physics D-Applied Physics* 41(7), 2008.
5. Huja, S.S., Hay, J.L., Rummel, A.M., and Beck, F.M., "Quasi-Static and Harmonic Indentation of Osteonal Bone," *Journal of Dental Biomechanics* 2010, 1–7, 2010.
6. Hay, J., "Complex Shear Modulus of Commercial Gelatin by Instrumented Indentation," Keysight Technologies, Inc., 2014, Document No: 5990-9745EN, Date Accessed: June 6, 2013; Available from: <http://cp.literature.keysight.com/litweb/pdf/5990-9745EN.pdf>.

## Nano Mechanical Systems from Keysight Technologies

Keysight Technologies, the premier measurement company, offers high-precision, modular nano-measurement solutions for research, industry, and education. Exceptional worldwide support is provided by experienced application scientists and technical service personnel. Keysight's leading-edge R&D laboratories ensure the continued, timely introduction and optimization of innovative, easy-to-use nanomechanical system technologies.

[www.keysight.com/find/nano](http://www.keysight.com/find/nano)

For more information on Keysight Technologies' products, applications or services, please contact your local Keysight office. The complete list is available at: [www.keysight.com/find/contactus](http://www.keysight.com/find/contactus)

### Americas

Canada	(877) 894 4414
Brazil	55 11 3351 7010
Mexico	001 800 254 2440
United States	(800) 829 4444

### Asia Pacific

Australia	1 800 629 485
China	800 810 0189
Hong Kong	800 938 693
India	1 800 112 929
Japan	0120 (421) 345
Korea	080 769 0800
Malaysia	1 800 888 848
Singapore	1 800 375 8100
Taiwan	0800 047 866
Other AP Countries	(65) 6375 8100

### Europe & Middle East

Austria	0800 001122
Belgium	0800 58580
Finland	0800 523252
France	0805 980333
Germany	0800 6270999
Ireland	1800 832700
Israel	1 809 343051
Italy	800 599100
Luxembourg	+32 800 58580
Netherlands	0800 0233200
Russia	8800 5009286
Spain	0800 000154
Sweden	0200 882255
Switzerland	0800 805353
	Opt. 1 (DE)
	Opt. 2 (FR)
	Opt. 3 (IT)
United Kingdom	0800 0260637

For other unlisted countries:  
[www.keysight.com/find/contactus](http://www.keysight.com/find/contactus)  
(BP-07-10-14)

This information is subject to change without notice.  
© Keysight Technologies, 2013 - 2014  
Published in USA, August 3, 2014  
5991-2630EN  
[www.keysight.com](http://www.keysight.com)